# CA05FL

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# **CA05FL**

## **Owner's Manual**

**Fuel Cell Technology  
Breath Alcohol Tester**

# Part Names


This device is a breath alcohol tester designed to measure the presence of alcohol in the body. Please make sure to read this manual thoroughly before use.

**DO NOT DRINK AND DRIVE**

# Display

Test result is displayed in BrAC / Promile data correspondent to the 3<sup>rd</sup> decimal readout at BAC

Measurement unit is set to the region standard.

"bAt" signal (  ) means low battery voltage.  
It is recommended to replace 1.5V Alkaline batteries.

This unit has Dual Color LCD Backlight.  
See the Error Notification section in the manual.

If the result is over 0.50‰ (or equivalent to it in other measurement standard) LCD backlight turns red to give notice.

# Operation

**The purchased product is already calibrated and can be used immediately.**

1. Insert a mouthpiece into the mouthpiece inlet.
2. Push the ON/OFF power switch under the display window.
  - The unit will show "Full display" along with a welcome sound. And then it starts countdown which will be taken about 15 ~ 35 seconds.
  - This warm-up process prepares the sensor and circuit for testing.
3. When you hear a beep sound & "bLo" display is on, please blow steadily into the mouthpiece until you would hear the "Beep" sound again. Stop blowing.
  - In this stage, if you don't blow within 20 seconds, it automatically turns off. Try to power on again.
4. After the circling line is on the display for 3 ~ 22 seconds, the test result (mg/l BrAC or ‰ Promile) will be displayed for 20 seconds along with beep sound.
  - If the test result is over 0.25mg/l / 0.50‰ (default setting) respectively, the LCD color will be changed to red along with "Alarm" sound.
  - The result will be displayed for 20 seconds. (light on for the first 15 seconds and light off for another 5 seconds). If you push the power switch during the last 5 seconds, the unit goes back to countdown stage for another test.
5. Finally, the unit will show "OFF" and turn off automatically along with a sound.
6. For the successive testing, try again from No.2 to No.4.
7. Space tests at least 2 minutes apart.

※ After using, you can turn the device off by pressing the power switch for 2 sec.

## Test Results Record (optional)

Latest 10 test results are stored in the memory.

To recall the test results, press the ON/OFF switch when "ALC" is on the display. (see above step 2.)

After showing "rEC", it will display the latest test result(d-1).

The next test results are displayed sequentially by pressing ON/OFF switch.

If the memory is full, it will erase the oldest data.


It will automatically go to countdown stage if you don't press the ON/OFF switch for 3sec.

# Error (Warning) Notification

Error Message	Cause & Action	Error Message	Cause & Action
	If the result is higher than 0.25mg/l BrAC (0.50‰ Promile), LCD Backlight turns red along with "ALARM" sound.		The battery voltage is too low to perform a test. The device will turn off automatically after showing this signal. Replace 1.5V Alkaline Batteries.
	If the result is above 2.00mg/l BrAC (4.00‰ Promile), "Hot" signal is displayed along with alarm sound.		The user had a short exhale time during the test, which likely means the user did not exhale steadily and continuously for at least 3~5 seconds. Push the power button to perform another test.
	When the test time is over 500, the unit shows "CAL". Even if the unit still works, it is required the re-calibration. Consult to the seller or local distributor for the service.		When the ambient temperature is out of working temperature range (5 °C ~ 40 °C) the unit shows "OUT" and then turn off. Keep and use the unit within the recommended temperature range.

\*\*\* N.B. : Any of these error occurs, the LCD Backlight turns red to give notice.


# Battery Replacement


※ Remove Battery tab before use

Apply a small amount of pressure on the battery cover located on the back of the unit and push(①) and away(②) from the unit to release the cover. Put AAA type Alkaline batteries as the following diagram.

Step 1


Step 2


※ Insert from negative direction


Step 3

Step 4

Step 5

\*\*\* N.B. : Make sure to insert the battery with the proper position of positive(+) and negative(-) as indicated.

# Precautions

1. Allow at least 20 minutes after eating or drinking before testing.
  - Any type of alcohol remaining in the mouth of saliva may interfere with proper testing and requires about 20 minutes to disperse. This includes alcohol-based mouthwash or similar products.
2. Avoid testing under strong wind or in the closed room with contaminated air.
3. When the "bAt" signal is on the display, it is advisable to replace 1.5V alkaline batteries. **DO NOT USE RECHARGEABLE BATTERIES.** When handling, pay close attention to reinserting the battery cap.
4. Do not blow cigarette smoke, food or liquid into the instrument because this may damage the sensor.
5. Do not use this unit as a tool to determine whether you or anyone else should operate any motor vehicle or machinery.
6. Keep the device out of reach of children!!!  
Keep it under 0 – 40 °C (32 – 104 °F). And this device is designed to be used in a temperature range of 5 – 40 °C (41 – 104 °F). Keep and use the unit within the recommended temperature range.
7. Correlation between breath alcohol content and blood alcohol concentration depends on many variables, including environmental factors (such as air quality, wind, humidity, temperature, etc.) and health conditions of the test subject.
8. In normal use, CA05FL requires re-calibration every 6 months or after 500 tests to obtain accurate results. Discuss with local distributors for re-calibration. If you do not let to calibrate the device, the device sensor may burn and device sensor warranty will void. This calibration is not included in the warranty, but is at the expense of the customer.
9. You cannot use the results of this unit in court of law. The driver holds full responsibility for driving under the influence of alcohol.
10. Do not use this device in conjunction with an ignition interlock device or while sitting in a motor vehicle .
11. An individual's testing result can continue to increase or decrease over time, depending on many factors. Obtaining a testing result reading at one time does not mean the testing result reading will be the same even a few minutes later.
12. Do not let yourself or anyone else drink and drive. Always use a designated driver when you or anyone else drinks alcohol.

# Specifications

Sensor	Advanced Fuel Cell Sensor
Accuracy	±0.005% BAC at 0.050% BAC
Detection Range	0.00 ~ 2.00mg/l BrAC 0.00 ~ 4.00‰ Promile
Warm-up Time	15 ~ 35 seconds
Response Time	3 ~ 22 seconds
Dimension	111mm x 35mm x 24mm
Weight	55g, (77g Battery included)
Power Supply	Two 1.5V "AAA" size Alkaline Batteries
Working Temperature	5 °C ~ 40 °C
Storage Temperature	0 °C ~ 40 °C
Calibration	Every 6 months or after 500 tests
Warranty	2 Year Warranty

## Warranty

Manufacturer warrants for 2 years the product to be free from defects in workmanship or material (excluding calibration) in normal service from the date of purchase.

Manufacturer's obligation under the warranty is limited the replacing, adjusting, or repairing the unit if returned along with proof of purchase.

This warranty is void if the unit has been tampered with or maliciously damaged.

Disposal: the product has to be disposed of in accordance with the Directive of the European Parliament and of the Council 2012/19/EU on disposal of electronic waste. The product is marked with a pictogram in accordance with EN 50419.

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